






FSH GENLISA™ ELISA

REF : KBD381

Ver 2.0


IVD

Enzyme Immunoassay for Quantitative Determination of FSH in serum and plasma.

| | | | |
|---|-----------------------------|---|--------------------------------|
| IVD | For In-vitro Diagnostic Use | REF | Catalog Number |
|  | Store At | LOT | Batch Code |
|  | Manufactured By |  | Biological Risk |
|  | Expiry Date |  | Consult Operating Instructions |

For In-vitro Diagnostic Use only. Purchase does not include or carry the right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of KRISHGEN Pudgala LLP is strictly prohibited.

REF KBD381

 96 tests

Krishgen Pudgala LLP Unit Nos#318/319, Shah & Nahar, Off Dr E Moses Road, Worli, Mumbai 400018. India. Tel: +91-22-49198700 | email: sales@krishgenpudgala.com

Introduction:

Follicle stimulating hormone is one of the gonadotrophic hormones, secreted by pituitary gland into the bloodstream. Follicle stimulating hormone is one of the hormones essential to pubertal development and the function of women's ovaries and men's testes. In women, this hormone stimulates the growth of ovarian follicles in the ovary before the release of an egg from one follicle at ovulation. It also increases estradiol production. In men, follicle stimulating hormone acts on the Sertoli cells of the testes to stimulate sperm production (spermatogenesis).

Intended Use:

The FSH GENLISA™ ELISA is intended for the quantitative determination of Follicle Stimulating Hormone in serum and plasma.

Principle:

The FSH GENLISA™ ELISA method employs sandwich enzyme linked immunosorbent assay (ELISA) technique. Monoclonal Anti-FSH Antibodies are pre-coated onto microwells. Standards/Samples/Controls are pipetted into microwells and Follicle Stimulating Hormone present in the sample is bound by the antibodies. Enzyme conjugate antibody is pipetted and incubated to form a complex. After washing microwells in order to remove any non-specific binding, the ready to use TMB substrate is added to microwells and color develops proportionally to the amount of Follicle Stimulating Hormone present in sample. Color development is then stopped by addition of stop solution. Absorbance is measured at 450 nm.

Materials Provided:

1. Microtiter Coated Plate (8x12 wells) - 1 no
2. Standard 0 (2 ml/vial) - 0 mIU/ml
3. Standards (0.5 ml/vial) - 5, 10, 20, 50, 100 mIU/ml
4. Control Serum - 0.5 ml
5. Enzyme Conjugate - 12 ml
6. (25X) Wash Buffer – 2 X 25 ml
7. TMB Substrate - 12 ml
8. Stop Solution - 12 ml
9. Instruction Manual.

Materials to be provided by the End-User:

1. Microtiter Plate Reader able to measure absorbance at 450 nm.
2. Adjustable pipettes and multichannel pipettor to measure volumes ranging from 25 ul to 1000 ul
3. Deionized (DI) water
4. Wash bottle or automated microplate washer
5. Graph paper or software for data analysis
6. Timer
7. Absorbent Paper

Handling/Storage:

1. Store main kit components at recommended storage temperature indicated on the component label.
2. Before using, bring all components to room temperature (18-25°C). Upon assay completion return all components to appropriate storage conditions.
3. The Substrate is light-sensitive and should be protected from direct sunlight or UV sources.

Health Hazard Warnings:

1. Reagents that contain preservatives may be harmful if ingested, inhaled or absorbed through the skin. Refer to the MSDS online for details.
2. To reduce the likelihood of blood-borne transmission of infectious agents, handle all serum and/or plasma in accordance with NCCLS regulations.

**Specimen Collection and Handling:**

Serum- Coagulate at room temperature for 10-20 minutes; centrifuge for 20-min at 2000-3000 rpm. Remove the supernatant. If precipitation appears, recentrifuge.

Plasma- Use EDTA or citrate plasma as an anticoagulant, mix for 10-20 minutes; centrifuge for 15-min at 2000-3000 rpm. Remove the supernatant carefully. If precipitation appears, recentrifuge.

Reagent Preparation:

1. Wash Buffer (1X) Dilution: To make Wash Buffer (1X), add 4 ml of Wash Buffer (25X) to 96 ml of DI water. This is the working solution.
2. Allow all components to reach RT (Room Temperature) prior to use in the assay.

Test Procedure:

1. All reagents should be allowed to reach room temperature before use.
2. Add **25 ul Samples, Standard, Control Serum** into appropriate wells.
3. Add **100 ul Enzyme Conjugate** in each well.
4. Incubate at room temperature for 90 minutes.
5. Aspirate and wash plate 5 times with **(1X) Wash Buffer** and blot residual buffer by firmly tapping plate upside down on absorbent paper. Wipe of any liquid from the bottom outside of the microtiter wells as any residue can interfere in the reading step. All the washes should be performed similarly.
6. Add **100 ul of TMB Substrate** to all wells.
7. Incubate at room temperature for 10-15 minutes.
8. Add **100 ul of Stop Solution**. Read result with an ELISA reader at 450 nm within 15 minutes of stopping the reaction.

Calculations:

Determine the Mean Absorbance (net of Blank) for each set of duplicate Standards and Samples. Using graph paper, plot the average value (absorbance 450 nm) of each standard on the Y-axis versus the corresponding concentration of the standards on the X-axis. Draw the best fit curve through the standard points. To determine the unknown FSH concentrations, find the unknown's Mean Absorbance value on the Y-axis and draw a horizontal line to the standard curve. At the point of intersection, draw a vertical line to the X-axis and read the FSH Concentration. If samples were diluted, multiply by the appropriate dilution factor.

Software which is able to generate a cubic spline curve-fit, 4-PL or a polynomial curve (2nd order) is best recommended for automated results.

Note:

It is recommended to repeat the assay at a different dilution factor in the following cases:

- If the sample absorbance value is below the first standard.
- If the sample reads more than 100 mIU/ml, then dilute with Standard 0. The obtained result should be multiplied by the dilution factor.

If the determination value is higher or lower than normal range, it means there is an abnormal result. The final result should be diagnosed in correlation with the clinical symptoms and other diagnostic methods.

Reference Value:

| Group | Phase | Range (mIU/ml) |
|---------|------------------|----------------|
| Males | - | 1.4-14 |
| Females | Follicular Stage | < 10 |
| | Midcycle Peak | 5-16 |
| | Luteal Phase | < 10 |
| | Postmenopausal | 25-150 |

It is recommended that each laboratory establishes its own normal and pathological reference ranges, as usually done for other diagnostic parameters, too. Therefore, the above mentioned reference values provide a guide only to values which might be expected.

Limitations of Method:

Any clinical diagnosis should not be based on the results of in vitro diagnostic methods alone. Physicians are supposed to consider all clinical and laboratory findings possible to state a diagnosis.

Performance Characteristics:**Sensitivity:**

Limit Of Detection: It is defined as the lowest detectable concentration corresponding to a signal of Mean of '0' standard plus 2*SD. 12 replicates of '0' standards were evaluated and the LOD was found to be 0.3 mIU/ml.

Specificity:

| Cross-reagent | Cross reactivity, % |
|---------------|---------------------|
| FSH | 100 |
| hCG | 0.001 |
| TSH | 0.0000 |
| LH | 0.019 |

Precision:

Intra-Assay-Variation

| Sample | N | Mean Value | Standard Deviation | CV (%) |
|--------|---|------------|--------------------|--------|
| 1 | 9 | 10.1 | 0.450 | 4.5 |

Inter-Assay-Variation

| Sample | Mean Value | Standard Deviation | CV (%) |
|--------|------------|--------------------|--------|
| 1 | 9.7 | 0.442 | 4.6 |

Safety Precautions:

- **This kit is For In-vitro Diagnostic Use only.** Follow the working instructions carefully.
- The expiration dates stated on the kit are to be observed. The same relates to the stability stated for reagents
- Do not use or mix reagents from different lots.
- Do not use reagents from other manufacturers.
- Avoid time shift during pipetting of reagents.
- All reagents should be kept in the original shipping container.
- Some of the reagents contain small amount of sodium azide (< 0.1 % w/w) as preservative. They must not be swallowed or allowed to come into contact with skin or mucosa.
- Source materials maybe derived from human body fluids or organs used in the preparation of this kit were tested and found negative for HBsAg and HIV as well as for HCV antibodies. However, no known test guarantees the absence of such viral agents. Therefore, handle all components and all patient samples as if potentially hazardous.
- Since the kit contains potentially hazardous materials, the following precautions should be observed
 - Do not smoke, eat or drink while handling kit material
 - Always use protective gloves
 - Never pipette material by mouth
 - Wipe up spills promptly, washing the affected surface thoroughly with a decontaminant.
- In any case GLP should be applied with all general and individual regulations to the use of this kit.



LIMITED WARRANTY

Krishgen Pudgala LLP does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the product; against defects in products or components not manufactured by Krishgen Pudgala LLP, or against damages resulting from such non-Krishgen Pudgala LLP made products or components. Krishgen Pudgala LLP passes on to customer the warranty it received (if any) from the maker thereof of such non-Krishgen made products or components. This warranty also does not apply to product to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by Krishgen Pudgala LLP.

THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of Krishgen Pudgala LLP shall be to repair or replace the defective product in the manner and for the period provided above. Krishgen Pudgala LLP shall not have any other obligation with respect to the products or any part thereof, whether based on contract, tort, strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall Krishgen Pudgala LLP be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of Krishgen Pudgala LLP with respect to the product. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

Krishgen Pudgala LLP, 2021.

THANK YOU FOR USING KRISHGEN PRODUCT!



Krishgen Pudgala LLP

Unit No.1/2, Om Sainath Commercial Complex,
Off Mankoli-Anjur Phata Road. Village Dapode, Bhiwandi 421302.












Regulatory Status:

| | |
|------------------|--------|
| CE Marked | Europe |
| FDA registered | USA |
| CDSCO registered | India |

SCHEMATIC ASSAY PROCEDURE

| | |
|---|---|
| 1 | All reagents should be allowed to reach room temperature before use. |
| 2 | Add 25 ul Sample, Standard, Control Serum into appropriate wells. |
| 3 | Add 100 ul Enzyme Conjugate in each well. Incubate at room temperature for 90 minutes. |
| 4 | Aspirate and wash plate 5 times with (1X) Wash Buffer and blot residual buffer by firmly tapping plate upside down on absorbent paper. Wipe of any liquid from the bottom outside of the microtiter wells as any residue can interfere in the reading step. All the washes should be performed similarly. |
| 7 | Add 100 ul of TMB Substrate to all wells. |
| 8 | Incubate at room temperature for 10-15 minutes . |
| 9 | Add 100 ul of Stop Solution . Read result with an ELISA reader at 450 nm within 15 minutes of stopping the reaction. |

SYMBOLS KEY

| | |
|---|-------------------------------|
|  | Microtiter Plate (8x12 wells) |
|  | Control Serum |
|  | Standards |
|  | Enzyme Conjugate |
|  | (25X) Wash Buffer |
|  | TMB Substrate |
|  | Stop Solution |
|  | Consult Instructions for Use |
|  | Catalog Number |
|  | Expiration Date |
|  | Storage Temperature |